



Kennebec
Pharmacy &
Home Care

CADD PRIZM VIP PUMP INFORMATION **FOR 5FU INFUSION PATIENTS**

- After approx. 1 minute, your screen will go blank to conserve the battery. You may bring it back at any time by pressing the “NEXT” key.
- When there is 5 mls remaining in your medication bag, the pump will beep and screen will alert you that “Reservoir Volume is Low”. Pump will keep running, but you will notice that both the green and the yellow light are now flashing. This is okay. Pump will perform this alert at each 1 ml dispensed after that. (i.e. at 4mls, 3mls, 2mls and 1ml). At the end of the infusion, the pump will give consecutive “double beeps” and screen will alert that “Reservoir Volume is Zero”. When this occurs, press “NEXT”, then remove the battery to turn pump off.
- If the pump beeps in the middle of the infusion, but appears to still be running. Check for the message “Battery Low” in the running screen. If this appears in the screen, change the battery, following the instructions “To Change Batteries During an Infusion” in your blue packet.
- If you have any questions about the pump or alarms, please call KPHC at 1-888-463-8083.