



CADD PRIZM VIP PUMP INSTRUCTIONS

INTERMITTENT MODE

TO ADD A NEW ANTIBIOTIC BAG

NOTE: Medication bag changes should be done when the infusion is in “KVO” mode, which is between doses -- not while the pump is giving a dose.

1. Stop pump by pressing the **STOP/START** key.

| |
|--------------------------------|
| Stop the Pump? Press Y or N |
|--------------------------------|

Press “Y”

2. Remove old battery and dispose of appropriately.
3. Remove old med bag tubing from pump.
4. Replace battery and battery compartment door. Pump will power up automatically and go through self test. (approx. 50 seconds.)
5. When self test is complete, pump will beep six times and come to the following screen:

PUMP DISPLAY

KEY PRESS

| |
|---|
| ***INT 6240D *** 9V Battery STOPPED Press NEXT to advance |
|---|

Attach the new tubing cassette. Pump will beep and register the following:

| |
|---|
| Admin set latched. Upstream sensor On. NEXT to continue |
|---|

Press “NEXT”

| |
|---|
| Reset Reservoir Volume to xxx.x ml? Press Y or N |
|---|

If this volume is correct to reset the total volume of your med bag, Press “Y”

| |
|-----------------------------------|
| Prime Tubing? Press Y or N |
|-----------------------------------|

Press “Y”

Disconnect tubing
from patient
Open clamps
Hold Y to prime

**Press and hold the “Y” to prime.
When primed, release the button.**

Continue Priming?

Press Y or N

Press “N”

Start the Pump?

Press Y or N

Press “Y”

Pump will start up and review the program parameters again, and display the following screen:

***INT 6240D ***
9V Battery KVO
Res Vol xxx.x ml
Press NEXT to advance

****green light will be blinking at the top of the pump while running and screen will go blank.**

PLEASE NOTE: Pump will begin to beep at 1ml increments beginning at 5mls remaining to notify of “Low Reservoir Volume”. This will be during a “KVO” phase and medication bag can be changed.