



## CADD PRIZM VIP PUMP INSTRUCTIONS

### *TPN MODE*

#### TO START TPN BAG

- **Remove battery door, insert 9Volt battery in correct position (pump will beep when battery is inserted correctly). Replace battery door.**
- **Pump will power up and automatically review program.**
- **Pump will beep six times when power up is complete.**
- **Attach tubing cassette to pump.**
- **Follow the displays below:**  
(\*\*\*\*Pump screen goes blank to conserve battery – press NEXT key to bring display back)

#### PUMP DISPLAY

#### ACTION

Reservoir latched.  
NEXT to continue

Press “NEXT”

Prime Tubing?  
Press Y or N

Press “Y”

Disconnect tubing  
from patient  
Open clamps  
Hold Y to prime

Press and hold “Y” key  
until fluid reaches the end of the  
tubing, then release.

Continue Priming?  
Press Y or N

If needing to prime more, press  
“Y”. If done, press “N”

Start the Pump?  
Press Y or N

Press “Y”

**\*\*Pump will go through the parameters again and start up. The screen will appear as follows:**

***TPN	6230C ***
9V Battery	RUNNING
Res Vol	xxxx.x ml
Press NEXT	to advance

**\*\*\* Green light will be blinking at the top of the pump and the screen will go blank.**

**When pump has completed the infusion or ramp down period (time-limited infusions only), it will beep several times, indicating that pump is entering the KVO (Keep Vein Open) mode. \*\*\*ONLY WHEN PUMP IS IN KVO, may you then proceed with the following to stop the pump:**

**PUMP DISPLAY**

**ACTION**

***TPN	6230C ***
9V Battery	RUNNING K
Res Vol	xxxx.x ml
Press NEXT	to advance

**Press "STOP/START"**

Stop the Pump?
Press Y or N

**Press "Y"**

Infusion Profile has been reset
------------------------------------

Reservoir Volume xxxx.x ml
<Resetting>

***TPN	6230C ***
9V Battery	STOPPED
Press NEXT	to advance

**Remove battery to turn pump off.**